

2019

UG 3rd Semester (Honours) Examination

PHYSICS

Paper - C6P

[Practical]

Full Marks : 20

Time : 3 Hours

*The figures in the margin indicate full marks.
Candidates are required to give their answers
in their own words as far as practicable.*

Distribution of Marks : Experiment - 15

LNB : 02

Viva : 03

Perform one experiment allotted through drawing cards.

1. Determine Mechanical Equivalent of heat by Callender and Barne's constant flow method.

(a) Working formula.

2

[Turn Over]

- (b) Circuit diagram. 1
- (c) Initial steady temperature difference of the thermometers. 1
- (d) Recording of steady temperatures at the inlet and outlet respectively, potential difference across the coil, current through it and rate of flow of liquid. 3×3

[for three different sets of currents and voltages]

- (e) Calculation of mechanical equivalent of heat. 1½
- (f) Accuracy. 1½

2. Determine thermal conductivity of copper by Searle's Apparatus.

- (a) Working formula. 2
- (b) Experimental diagram. 1
- (c) Initial steady temperatures recorded by thermometers placed at the inlet and outlet of coiled tube; and two other points on the bar. 1
- (d) Recording of temperatures of four thermometers at steady state during steam flow. 4
- (e) Collection of water in given time 't'. 2

- (f) Diameter of the bar by slide callipers. (V. C. supplied) 2
- (g) Distance between the thermometers on the bar. 1
- (h) Calculation 1
- (i) Accuracy 1
3. Determine thermal conductivity (k) of copper by Angstrom's method.

Apply a periodic square heat pulse to one end of a copper rod while other end remains in room temperature. Record temperatures by two thermistors or other devices at a small distance Δx apart, as a function of time. Apply Fourier transformation for two distinct measurements of the thermal conductivity of the copper rod for 1st and 3rd harmonics of heat wave.

- (a) Theory 2
- (b) Arrangement for square heat pulse. 1
- (c) Time-temperature data collection with proper acquisition unit. 5

[Turn Over]

(d) Fourier analysis of data 5

(e) K for 1st and 3rd harmonics of heat wave. 2

4. Determine thermal conductivity of a bad conductor by Lee and Charlton's disc method.

[Mass and thickness of the lower disc, Diameter and thickness of the bad conducting disc and steady state temperature are to be supplied.]

(a) Working formula with Bedford's correction 2

(b) Time-temperature record during cooling (after direct heating by steam) 4

(c) Time-temperature graph to find rate of cooling 3

(d) Bedford's correction 1

(e) Calculation 2

(f) Accuracy 1

(g) Mention the process of measuring thickness of the experimental disc and the process of recording steady state temperature. 1+1

5. Determine the Temperature coefficient of Resistance (α) by Platinum Resistance Thermometer (PRT)

[Measure resistance of PRT at two different known temperatures preferably at ice and steam]

- (a) Working formula and circuit diagram 1+1
- (b) Data for electrical mid point. 1
- (c) Data for determination of the resistance of the PRT ice and steam 4+4
- (d) Evaluation of resistance per unit length (ρ) of the bridge-wire. 1
- (e) Calculation of resistance at two different temperatures. 1+1
- (f) Evaluation of temperature coefficient (α) of resistance. 1

6. Study the variation of Thermo-Emf of a Thermocouple with difference of temperature of its two Junctions.

[Resistance of the potentiometer wire is to be supplied]

[Turn Over]

- (a) Working formula and circuit diagram 1+1
- (b) Calculation of R to be put in series with potentiometer wire for $5 \mu\text{V}/\text{cm}$ drop 1
- (c) Datas for $e-t$ graph (at least six points) + calculation of e 6+2
- (d) $e-t$ graph 2
- (e) Determination of thermoelection power 1
- (f) Accuracy 1

7. Calibration of a thermocouple within the temperature range of 80°C to 40°C with cold junction at ice bath. Hence finding melting point of Wax using calibration curve (null point length (l) vs. temperature (t) graph). (Rp is to be supplied]

- (a) Working formula t circuit diagram 1+1
- (b) Calculation of R to be put in series with potentiometer wire fo r $5 \mu\text{V}/\text{cm}$ drap, considering given Rp. 1
- (c) Data for null-point length (l) vs. temperature graph 3

- (d) Drawing of null point length (l) vs. temperature (t) graph. 2
- (e) Data for null point length vs. time graph during melting or freezing of wax. 3
- (f) Null point length vs. time graph. 2
- (g) Determination of melting point using calibration curve. 1
- (h) Accuracy 1
8. Calibration of a thermocouple within the temperature range of 80°C to 40°C with cold junction at ice bath, by direct measurement using OPAMP.
- (a) Theory and circuit diagram 1+1
- (b) OFF-SET null adjustment 2
- (c) Hot junction temperature vs output voltage data for at least six difference temperatures. 6
- (d) Calibration curve (hot junction temperature vs. output voltage) 2
- (e) Calculation of thermo-electric power from the calibration graph 2
- (f) Accuracy 1